Technical Specifications



- Overall dimensions
- Mattress platform
- Hi-Low
- Side rail height
- Side rail length
- Backrest angle
- Kneerest angle
- Bed weight
- Safe working load
- Casters

2260 mm × 1010 mm 2000 mm × 810 mm 770 mm - 480 mm 410 mm 1620 mm 0 - 70° 0 - 35° 135 kg

220 kg

125 mm

Technical and Quality Assurance Standards IEC 60601-2-52 ISO 13485 ISO 14971





Positions







B3-U5-7, Solaris Dutamas, Jalan Dutamas 1, 50480 Kuala Lumpur, Malaysia Tel: (+60 3) 62113120, 62112120, 62112120 www.piyatech.com email: info@piyatech.com



General Bed 35000MC1 Type 2

MDA Establishment License No.: MDA-1216-Kl20 MDA Product Reg. No.: GMD45614744918A



Features



Side Rails

The parallel side rails of the 35000MC1 offer a unique and stylish design produced from a combination of blow molded plastic and stee I bars to provide a protective height and length. The side rails adhere to the highest standards of safety to prevent patient falls and entrapment.



Side Rails

the side rails to fold down.

Locking Mechanism Pushing up the red lever on either side of the bed engages the release mechanism allowing



Bed Ends

The seamless blow-molded bed ends are lightweight, durable, and easy to clean. The bed ends are of high quality and compatible with steam cleaning. To have better access to the patient's head in emergency cases, the caregiver can easily remove the bed ends by grasping the hand holds and lifting them up in one single step. Moreover, the bed ends have integrated steel fittings which make them extremely strong so the nursing staff can move the bed by holding and pushing them without worries.



Mattress Platforms

The mattress platforms exhibit exceptional qualities, including smooth surfaces, high durability, and effortless removal for easy cleaning and disinfection. Additionally, the mattress platforms are designed to be compatible with steam cleaning, ensuring thorough and effective sanitation.



Foldable Cranks

Back, knee, and hi-low positions can be adjusted by three foldable crankslocated at foot end of the bed.





Urine Bag Holder Beds are equipped with a urine bag holder which can be easily installed on either side of the bed.

General Bed 35000MC1 Type 1





Casters

125 mm casters run extremely smoothly; the quality that makes it significantly easy for the nursing staff to transport the bed in the medical centers. Having diagonal locks, the casters can be locked to prevent the bed from unintended movement.



Bumpers

Innovative rotary bumpers which are located on four corners of beds can prevent to a large extent the damages that may cause collision of the bed to the surroundings during transportation.

